

पेटेंट कार्यालय
शासकीय जर्नल

**OFFICIAL JOURNAL
OF
THE PATENT OFFICE**

निर्गमन सं. 27/2023
ISSUE NO. 27/2023

शुक्रवार
FRIDAY

दिनांक: 07/07/2023
DATE: 07/07/2023

पेटेंट कार्यालय का एक प्रकाशन
PUBLICATION OF THE PATENT OFFICE

(12) PATENT APPLICATION PUBLICATION

(21) Application No.202311038333 A

(19) INDIA

(22) Date of filing of Application :05/06/2023

(43) Publication Date : 07/07/2023

(54) Title of the invention : DISPOSABLE HAEMORRHOID DEVICE FOR THE TREATMENT OF HAEMORRHOIDS

(51) International classification :A61B 171200, A61K 090000, A61K 473400, A61P 031000, A61P 091400
(86) International Application No :NA
Filing Date :NA
(87) International Publication No : NA
(61) Patent of Addition to Application Number :NA
Filing Date :NA
(62) Divisional to Application Number :NA
Filing Date :NA

(71)Name of Applicant :
1)IQBAL HUSEN SHAIKH
Address of Applicant :IQBAL HUSEN SHAIKH SEGURO
LIFESCIENCES PLOT 2(i)A EPIP PHASE 1 JHARMAJARI ----

Name of Applicant : NA
Address of Applicant : NA
(72)Name of Inventor :
1)IQBAL HUSEN SHAIKH
Address of Applicant :IQBAL HUSEN SHAIKH SEGURO
LIFESCIENCES PLOT 2(i)A EPIP PHASE 1 JHARMAJARI ----

(57) Abstract :

ABSTRACT DISPOSABLE HEMORRHOID TREATMENT DEVICE This invention is directed to a medical device. A Disposable Haemorrhoid Device is the equipment used to achieve low temperature <35° C for treatment of Hermorrhoids. Device is designed for the direct application of low temperature therapy. It is designed with consideration of human anatomy and medical concepts for treatment of this medical condition. The Device has a first section which is smooth blunt tip in bullet shape, a middle section which is in cylindrical body and end section is in oval disk shape with a flap at right angle for holding the device and apply pressure to insert the first section into the anal orifice. A Non Toxic, low thermal conductive Gel is filled in the device. For usage The Haemorrhoid Device is kept in a refrigerator Freezer at temperature -15° C to -18° C for 90 min. before use. This device is for onetime use. By this device the direct application of low temperature provides immediate relief to pain, itching and inflammation and Haemorrhoids.

No. of Pages : 8 No. of Claims : 8